AxonGram

Reagenzien zur Gramfärbung
Réactifs selon la coloration de gram
Reagents according to the coloration of gram

Packungsprospekt
Prospectus d'emballage
Package insert

R1	500 ml	REF	AL0100
R2	500 ml	REF	AL0101
R3	500 ml	REF	AL0102
R4	500 ml	REF	AL0103

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Intended use

The reagents are intended for Gram's staining of bacteria and fungi.

Test principle

The crystal violet oxalate **R1** stains the cell walls violet. The crystal violet is fixed with Lugol's solution **R2** and the gram-negative germs are decolorized with the acetone ethanol mixture **R3**. The fuchsin solution **R4** serves as the counterstaining agent (gram-negative germs appear red and gram-positive germs remain violet).

Reagents/Composition

Crystal violet solution

REF AL0100	Flammable Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment	500 ml	Crystal violet 2% Ammonium oxalate 0.8% Ethanol 16% Demineralized water ad 500 ml
---------------	--	--------	--

Lugol's solution

REF AL0101	R2	500 ml	Iodine 1.3% Potassium iodide 2% Polyvinylpyrrolidone 5% Demineralized water ad 500 ml
---------------	-----------	--------	--

Decolorizing solution

REF AL0102	R3  F Highly flammable Keep away from sources of ignition – No smoking  Xi Irritant Irritating to eyes	500 ml	Acetone 25% Ethanol 75%
---------------	---	--------	----------------------------

Fuchsin solution

REF AL0103	R4	500 ml	Basic Fuchsin 0.1% Demineralized water ad 500 ml
---------------	-----------	--------	---

Stability

At 15 to 25°C, the reagents are stable until the expiration date mentioned on the labels, even after opening of the vials.

Analysis material/Sample material

Both bacteriological samples and cultivated micro-organisms can be used.

How to perform the test

1. Apply the analysis material on a clean slide.
2. Let it dry.
3. Fix with heat.
4. Staining:
 - Cover the preparation with a layer of crystal violet solution **R1** for one minute.
 - Carefully rinse with water.
 - Cover with a layer of Lugol solution **R2**, let drip and cover again with Lugol solution **R2**, then let absorb for one minute.
 - Carefully rinse with water and decolorize thoroughly with decolorizing solution **R3**.
 - Carefully rinse with water and cover with a layer of fuchsin solution **R4** for one minute. Then carefully rinse with water and let dry.
5. Reading and interpretation under the microscope (oil immersion with 1,000x magnification).

Quality control

Quality control is recommended with the following reference strains:

Escherichia coli	ATCC 25922	Gram-negative rods: red coloring
Staphylococcus aureus	ATCC 25923	Gram-positive cocci: violet coloring

Limitations/Interferences

The Gram's stain should be performed with fresh cultures (24 to 48 h). Older cultures decolorize more easily. In case of antibiotic therapies, gram-positive germs decolorize very rapidly.

Performance characteristics

In an external study, performance data were ascertained with 150 gram-positive and gram-negative germs, yeasts and molds. All gram-positive and gram-negative germs were precisely and thoroughly identified.

Waste disposal

It is the laboratory's responsibility to comply with the respective national regulatory requirements when disposing of the analysis material and the reagents (see Material Safety Data Sheets).

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Axon Lab AG