# CountAX TPC / RB

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# CountAX TPC / RB

## Intended use


Germ indicator for determining the total plate count (TPC) and detecting yeasts and mold (RB) on solid and in liquid materials. The germ indicator can also be used for transporting samples.

## Package contents

- 10 dip slides
- 10 labels

## Basic principle

The CountAX TPC/RB Germ Indicator comprises 2 culture media:

Side 1 with the TPC agar can be used to rapidly grow the most common bacteria and fungi. Side 2 with the RB agar can be used to grow yeasts and mold. Antibiotics are added to inhibit bacteria growth.

## Composition

Total plate count agar (TPC), yellowish	Rose-Bengal agar (RB), magenta
Casein peptone	Casein peptone
Soy peptone	Yeast extract
Sodium chloride	Glucose
Lecithin	Bengal red
Polysorbate 80	Antibiotics
Sodium thiosulfate	Disodium phosphate
L-Histidine	Magnesium sulfate
Agar	Agar
pH 7.3	pH 7.2

## Storage

The CountAX is stored at 15–25°C. Temperature fluctuations and drafts must be avoided. Do not freeze. Do not use the product after the expiration date (EXP) specified on the package.

Do not use the product if it becomes discolored or dried out, if the agar layers become detached from the plastic medium, or if there are indications of bacterial/fungal growth. Do not touch the growing colonies as they might contain pathogenic germs.


## Sampling

Any contamination of the agar surfaces (growing media) must be avoided. The agar surfaces must be placed with full contact on the analysis material. Once the sample has been taken, screw the tube tightly shut again.

## Performing the test


1. Unscrew the CountAX tube and remove the culture medium without touching the agar surfaces.

2. Contact inoculation

The flexible holder with the bracket enables easy handling (a). Solid surfaces can be tested by pressing each side of the germ indicator firmly against the relevant surface and holding it steadily for 3–5 seconds (b).

2a. Dipping process

Liquid samples are tested by dipping the germ indicator into the liquid for three to four seconds until the agar surfaces are completely covered (c). Use filter paper to wipe off drops from the sample.

2b. Smearing process

Semi-solid materials or objects that are difficult to access can be tested by smearing a sterile swab carefully over an area that has an edge surrounding it, for example. If the object is dry the swab should first be moistened with sterile water. A moistened swab can also be used to obtain samples from powders (e.g. spices) or viscous fluids.

After smearing the sample area, roll the swab carefully across the germ indicator's agar surfaces from left to right and from bottom to top (d).

3. Insert the culture medium again in the tube and screw the tube tightly shut.

4. Complete the sample identification label and affix it to the tube.

5. Incubate the germ indicator tightly sealed in its tube for 3–5 days at 22–25°C.

## Transporting the inoculated CountAX tubes

The inoculated germ indicator can be transported at 15–25°C.

## Interpretation of the results

After incubation, remove the germ indicator from its tube, calculate the germ count (number of colony-forming units, CFU) and examine the color reactions.


Yeasts and mold grow on the Rose Bengal medium. In this case, certain bacteria are inhibited because antibiotics have been added.

The following threshold values may be considered as a basis for evaluating the overall contamination level in food:

Good	$< 10^5$ CFU/g
Permissible	$10^5 - 10^7$ CFU/g
Poor	$> 10^7$ CFU/g

The threshold values for cooked foods are lower than those for raw foods.

The presence of yeasts and mold in foods always indicates that the product has been handled wrongly or a lack of hygiene.

It is difficult to specify threshold values for the total germ count on work and other surfaces as the criteria depend on usage. However, a surface may be regarded as being inadequately cleaned if it has 5 colonies/cm<sup>2</sup> after cleaning. This figure corresponds roughly to 45 colonies per slide on the TPC agar side.

The following threshold values may be considered as a rough basis for evaluating the level of contamination:

Contact inoculation		
Clean	0	CFU/side
Contaminated	1–10	CFU/side
Highly contaminated	> 10	CFU/side

### Limitation of method

When using the CountAX germ indicator as a contact slide, the system is identical to the contact plate procedure in terms of sensitivity, whereas the dipping and smearing processes have a detection limit of 1,000 CFU/ml. The permissible TPC for drinking water is too low to be detected using CountAX slides.

Only results that have been achieved using the same inoculation systems should be compared. Keep away from light. Store the dip slides in the dark to prevent photo-oxidation of the dye.

### Quality control

Growth characteristics of the medium can be tested with the following strains:

	TPC	RB
<i>Candida albicans</i> ATCC 10231	+	+
<i>Aspergillus niger</i> ATCC 16404	+	+
<i>Escherichia coli</i> ATCC 25922	+	–

(+) Growth (–) Partial or complete inhibition

### Waste disposal

All laboratories are responsible for proper disposal of waste in accordance with the legal requirements.

### References

1. J.C.G. Ottow: *Mycologia* 64, p. 304–315 (1972)
2. F.P. Downes and Ito K. (Eds.): 2001, *Compendium of Methods for the Microbiological Examination of Foods*, 4th Ed. APHA, Washington D.C.

**Abkürzungen/Abréviations/Abbreviazioni/Abbreviations  
Symbole/Symboles/Simboli/Symbols**

Symbol/Symbole/Simbolo	Bedeutung/Signification/ Significato
	Bestellnummer / Numéro d'article / Numero di catalogo / Catalogue number
	Packungsinhalt / Contenu du coffret / Contenuto della confezione / Contents of kit
	Hersteller / Fabricant / Fabbricante / Manufacturer
	Lagertemperatur / Température de conservation / Temperatura di conservazione / Storage temperature
	Verwendbar bis / Date de péremption / Data di scadenza / Use before
	Chargenbezeichnung / Numéro du lot / Codice del lotto / Lot number
 <a href="http://www.axonlab.com">www.axonlab.com</a>	Gebrauchsanweisung beachten / Consulter le mode d'emploi / Leggere le istruzioni per l'uso / Consult the instructions for use
	Einmalgebrauch / Usage unique / Per uso unico / single use
	Vor Licht schützen / français / italienisch / english

### **Schweiz**

Axon Lab AG  
Täferstrasse 15  
CH-5405 Baden  
Tel. +41 (0)56 484 80 80  
Fax +41 (0)56 484 80 99  
info@axonlab.ch  
www.axonlab.ch

### **Suisse**

Axon Lab SA  
En Budron E9  
CH-1052 Le Mont-sur-Lausanne  
Tél. +41 (0)21 651 71 50  
Fax +41 (0)21 651 71 59  
info\_f@axonlab.ch  
www.axonlab.ch

### **Deutschland**

Axon Lab AG  
Heinrich-Otto-Strasse 1  
D-73262 Reichenbach / Stuttgart  
Tel. +49 (0)7153 9226-0  
Fax +49 (0)7153 9226-26  
info@axonlab.de  
www.axonlab.de

### **Österreich**

Axon Lab AG  
Gewerbepark 1  
A-6404 Polling  
Tel. +43 (0)523 88 77 66  
Fax +43 (0)523 88 77 66-12  
info@axonlab.at  
www.axonlab.at