UrinAX CL/MC/NOR

Transportmedium, Keimzahlbestimmung und selektive Anzucht von Bakterien im Urin
Milieu de transport, dosage des germes et développement sélectif de bactéries dans l'urine
Mezzo di trasporto, conta batterica e moltiplicazione selettiva dei batteri nell'urina
Transport medium, bacterial count and selective propagation of bacteria in urine

Gebrauchsanweisung Notice d'utilisation Istruzioni per l'uso Instructions for use

REF AL0204 Cont.

- 10 Dip-Slides / Lames immergées / Piastra a immersione / Dip slides
- 10 Patientenetiketten / Etiquettes pour patient / Etichette paziente / Patient labels

In vitro diagnosticum
Usage in vitro
Per uso diagnostico in vitro
For in vitro diagnostic use

Axon Lab AG, CH-5405 Baden
Tel. +41 (56) 484 80 80
Fax +41 (56) 484 80 99
info@axonlab.ch · www.axonlab.ch

Schweiz
Axon Lab AG
Täferstrasse 15
CH-5405 Baden
Tel. +41 (0)56 484 80 80
Fax +41 (0)56 484 80 99
info@axonlab.ch
www.axonlab.ch

Suisse
Axon Lab SA
En Budron E9
CH-1052 Le Mont-sur-Lausanne
Tél. +41 (0)21 651 71 50
Fax +41 (0)21 651 71 59
info_f@axonlab.ch
www.axonlab.ch

Deutschland
Axon Lab AG
Heinrich-Otto-Strasse 1
D-73262 Reichenbach/Stuttgart
Tel. +49 (0)7153 9226-0
Fax +49 (0)7153 9226-26
info@axonlab.de
www.axonlab.de

Österreich
Axon Lab AG
Gewerbepark 1
A-6404 Polling
Tel. +43 (0)523 88 77 66
Fax +43 (0)523 88 77 66-12
info@axonlab.at
www.axonlab.at

UrinAX CL/MC/NOR

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6

Field of application

Transport medium, bacterial count and selective propagation of bacteria in urine.

Package contents

- 10 UrinAX dip inoculum culture media
- 10 patient labels

Principle

UrinAX consists of a sealed container with 3 culture media:

The CLED culture medium is used to determine the total bacterial count in the urine. The MacConkey Agar is used for the selective propagation of Enterobacteriaceae. The agar with norfloxacina (NOR) is used for testing sensitivity to this substance.

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

10

Interpretation of results

After incubation, assess the bacteria growth on CLED Agar. Compare the bacterial count on the CLED Agar side with the following readout:

Bacterial count on CLED Agar/ml

Notes: Bacterial counts 10^4 bacteria/ml are usually not considered pathological. Bacterial counts between 10^4 and 10^5 are an uncertain result and should be repeated. Bacterial counts of >math>10^5</math> bacteria/ml indicate a urinary tract infection.

A mixture of different bacteria strains on the UrinAX is very probably due to contamination of the urinary specimen.

MacConkey Agar:

Assess the presence of colonies. With a positive lactose reaction, the colonies are red. With a negative lactose reaction, the colonies are colorless. Carry out biochemical or immunological identification of the bacteria.

NOR Agar:

Assess the presence of bacteria growth. Bacteria growth on this agar indicates a resistance to norfloxacina (NOR).

Quality control

Growth characteristics of the medium can be tested with the following strains:

	CL	MC	NOR
Escherichia coli ATCC 25922	+	+	-
Staphylococcus aureus ATCC 25923	+	-	-
Proteus mirabilis ATCC 12453	+	+	-
Enterococcus faecalis ATCC 19433	+	-	+

(+) Growth (-) Partial or complete inhibition

Disposal

All laboratories are responsible for proper disposal of waste in accordance with the legal requirements.

References

1. ASM Approved Guidelines – 1996: Guidelines for Assuring of Medical Microbiological Culture Media
2. Société Française de Microbiologie: Recommendations 2007
3. Westenfelder, M. et al. : Infection 15, 20-24 (1987)

9